

Murisphere Step-by-Step Tutorial

Murisphere Step-by-Step Tutorial

Audience: Technicians, PIs, and Facility Admins

1. Sign In

Open Murisphere in a browser and log in with your assigned credentials.

Technician: routine room operations

PI: project and forecasting oversight

Admin: configuration, import, and governance

2. Set Phone-Reachable Scan URL

Go to

Cages

tab and set

Scan Base URL

.

Example:

http://10.0.0.24:8000

or your production domain.

3. Generate and Print Cage Cards

Click

Generate + Print

to produce cage cards with barcode links.

Allow browser pop-ups so the print dialog opens.

4. Scan in the Animal Room

How to use a phone to scan a printed cage card:

Hold the printed card flat under room lighting.

Open your phone's default Camera app (iPhone Camera or Android Camera).

Point the camera at the printed barcode until a link banner appears.

Tap the banner/link to open Murisphere in your mobile browser.

Confirm the cage ID on screen matches the printed card.

Tap

Open Murisphere (login for edits)

to move into authenticated edit mode.

If scanning fails:

move closer, reduce glare, and verify the printed URL uses your phone-reachable host (not localhost

).

5. Perform Cage Updates

In

Scan/Edit

tab, verify cage identity and update:

Male/Female counts

Breeding status

Notes

Click

Save Changes

. Updates are logged in audit history.

6. Record Lifecycle Events

Use lifecycle and breeding forms for:

Weaning

Transfers

Timed mating and plug checks

(truncated; see HTML tutorial for full detail)