

Murisphere Step-by-Step Tutorial

Murisphere Step-by-Step Tutorial

Audience: Technicians, PIs, and Facility Admins

1. Sign In

Open Murisphere in a browser and log in with your assigned credentials.

Technician: routine room operations

PI: project and forecasting oversight

Admin: configuration, import, and governance

2. Set Phone-Reachable Scan URL

Go to

Cages

tab and set

Scan Base URL

.

Example:

http://10.0.0.24:8000

or your production domain.

3. Generate and Print Cage Cards

Click

Generate + Print

to produce cage cards with barcode links.

Allow browser pop-ups so the print dialog opens.

4. Scan in the Animal Room

Use your phone camera to scan the barcode on a cage card.

The scan opens cage info in browser immediately. Tap

Open Murisphere (login for edits)

to continue with authenticated editing.

5. Perform Cage Updates

In

Scan/Edit

tab, verify cage identity and update:

Male/Female counts

Breeding status

Notes

Click

Save Changes

. Updates are logged in audit history.

6. Record Lifecycle Events

Use lifecycle and breeding forms for:

Weaning

Transfers

Timed mating and plug checks

Harvest and retirement schedules

7. Use Reports and Imports

In

Reports

tab:

Export CSV/XLSX/PDF

Upload genotyping CSV

Import onboarding Excel (Admin)

8. Monitor Compliance and Capacity

(truncated; see HTML tutorial for full detail)