

Murisphere Role-Based Tutorial

Technician, Facility Manager, and Researcher Workflows Across Workspace and Chat

Murisphere Operations

2026-04-13

Contents

Purpose	3
Training Setup	3
Create the tutorial database	3
Training accounts	3
What this seed contains	3
Real examples used in this guide	4
How To Use This Guide	4
Best device mix	4
Which mode should you use?	4
Shared rule for all roles	4
Biological Background You Learn Along The Way	4
Why timing matters	4
Why genotype matters	5
Why pedigree matters	5
Why compliance fields matter biologically	5
Shared Foundations For All Roles	5
What a cage card must do	5
Phone QR workflow	6
Common prompts used across all roles	6
Role 1: Technician	7
What a technician cares about daily	7
Typical technician data entries	7
Typical technician reports	7
Technician workflow in the workspace	8
1. Start from the dashboard or Cages	8
2. Generate a cage card	8
3. Scan in the room and reopen the cage	8
4. Use the workspace for visual review	10
Technician workflow in chat	10
1. Start the shift with the morning brief	10
2. Open the same cage directly	10
3. Notice the abnormal condition before editing	10

4. Attempt the update and observe the hard-stop	10
5. Still do useful work: complete the due task	10
6. Re-check the queue	11
7. Ask for the weaning queue	11
8. Ask for mortality follow-up	11
9. Batch-print the next room pass	12
Technician overlap you should not skip	12
Technician mini-mission	12
Role 2: Facility Manager	12
What a facility manager cares about daily	12
Typical manager data entries	12
Typical manager reports	13
Facility manager workflow in the workspace	13
1. Start from analytics and compliance	13
2. Move into analytics for capacity and flow	14
3. Generate exports	14
Facility manager workflow in chat	14
1. Start with the morning brief	14
2. Ask the next obvious question	14
3. Check room pressure	15
4. Pull chargeback or SLA context	15
Facility manager overlap you should not skip	16
Facility manager mini-mission	16
Role 3: Researcher / PI	16
What a researcher cares about daily	16
Typical researcher data entries	16
Typical researcher reports	17
Researcher workflow in the workspace	17
1. Use the workspace for samples, cohorts, and planning	17
2. Use pedigree when the biology matters	17
3. Use planner views when demand is changing	18
Researcher workflow in chat	18
1. Open the main project	18
2. Move from project to a specific cage	18
3. Ask about reports or readiness	18
4. Use the sample/result example	19
Researcher overlap you should not skip	20
Researcher mini-mission	20
Shared Cross-Role Workflows	20
Workflow: Print a cage card	20
Workspace path	20
Chat path	20
Workflow: Scan a printed card with a phone	20
Workflow: Move between modes without losing context	20
Fun Applications To Learn Together	21
1. Weekly colony huddle	21
2. Facility operations check-in	21

3. Research planning conversation	21
4. Pedigree teaching moment	21
Troubleshooting	21
Final Checklist	21

Purpose

This tutorial teaches Murisphere through the three roles that actually drive mouse colony work every day:

1. **Technician**
2. **Facility manager / vivarium administrator**
3. **Researcher / PI**

Murisphere now supports **two full operating modes**: - a **traditional workspace** for visual review, dashboards, lists, reports, and batch actions - a **chat-first console** for direct, phone-friendly commands and quick operational work

The important point is not to choose one forever. The important point is to use the mode that is fastest and safest for the task in front of you.

This tutorial is built from a **real seeded training dataset** generated on **April 13, 2026**. The cage codes, projects, sample records, planner scenarios, and alerts below come from that training dataset.

Training Setup

Create the tutorial database

```
python3 -m venv .venv
source .venv/bin/activate
pip install -r requirements.txt
./venv/bin/python seed_tutorial_demo.py --db training_demo.db --force
MURISPHERE_DB=training_demo.db ./venv/bin/python app.py
```

Training accounts

- Admin / facility manager: `admin@murisphere.local` / `admin1234`
- Technician: `tech@murisphere.local` / `tech1234`
- Researcher / PI: `pi@murisphere.local` / `pi1234`

What this seed contains

Category	Count
Labs	20
Cages	3,000
Projects	73
Animals	372
Litters	48
Breeding pairs	48
Sample records	48
Planner scenarios	6
Alert-driving tasks	80

Category	Count
----------	-------

Real examples used in this guide

Topic	Example
Technician breeding cage	F1-L01-C0006
Technician larger breeder cage	F1-L01-C0008
Researcher sample/result cage	F1-L01-C0012
Resulted sample	SMP-0004 on F1-L01-C0012-P01
Main project example	L01-PRJ-01
Main planner example	Neurogenetics Lab Cohort Plan
Non-expired protocol example lab	Developmental Signaling Unit

How To Use This Guide

Best device mix

The ideal training setup uses: - a **phone** for scanning printed cage cards - a **laptop or tablet** for the workspace - an optional second browser tab for the **chat console** at `/chat/`

Which mode should you use?

Mode	Best for	Typical examples
Workspace	dashboards, list views, analytics, reports, batch printing, visual alerts	room utilization, compliance review, project inspector, genotyping workspace
Chat	direct actions, quick lookup, phone work, intent-driven updates	Open cage F1-L01-C0006, Show overdue tasks, Print cage card for F1-L01-C0006
Mixed	moving fluidly between both	print in workspace, scan with phone, update in chat, return to workspace for review

Shared rule for all roles

Print the cage card -> scan the QR with a phone -> open the cage in the browser -> choose the fastest safe mode -> complete the task immediately.

Biological Background You Learn Along The Way

Murisphere is not just a logistics tool. The fields on the card and in the app are meaningful because they affect breeding, welfare, and experimental interpretation.

Why timing matters

- **DOB** drives weaning windows, age-matched cohorts, and breeding readiness.
- **Date of weaning (DoW)** matters because delayed or early weaning changes both welfare and downstream study timing.

- **Plug checks, pair setup, and breeder age** change whether a line is actually productive.

Why genotype matters

A genotype is not only a label. It is an inheritance and experiment-readiness question.

For example: - a **Cre/+** animal may be useful for one project but not another - a **f1/f1** animal may need to be paired with the right driver line - missing genotype results can stall a project even when the colony appears numerically healthy

Why pedigree matters

If a line behaves unexpectedly, the pedigree helps answer: - who the sire and dam were - whether the litter pattern is plausible - whether poor survival or odd genotype mix is isolated or repeated

Why compliance fields matter biologically

Protocol, welfare, veterinary, and quarantine fields are not only administrative: - they change what can be done to the cage - they can stop in-room updates or project handoffs - they protect both animal welfare and study validity




Figure 1: Pedigree and lineage view used to explain biological context

Shared Foundations For All Roles

What a cage card must do

A cage card is not decoration. It is the physical hand-off between room work and the database.

A complete Murisphere card should show: - cage code - PI / lab - project codes - room / rack / slot - protocol number, description, and expiration date - breeding status - strain and genotype summary - male

/ female / total counts - tracked animals - litter history - QR code that opens the cage in the browser

Murisphere Cage Card (Rendered Preview)

Designed for fast room-floor scanning, complete husbandry context, and compliance readiness

Cage: F1-L03-C0215
 Protocol: IACUC-2026-0312

Group Owner: Dr. Meosha Hudson
 Description: PdxCre longitudinal cohort

Group Name: Hudson Lab
 Breeding: Pairing 35 (Timed Mating)

Projects: L03-PRJ-01, L03-PRJ-05
 Protocol Expires: 12/31/2026

Location: Room 1, Breeding Rack 1, Slot 6A
 Population: M1 / F2 / T3

Tracked IDs shown: 3 of 3



Scan with phone camera

ID	Sex	DOB	Genotype	Status
258	M	10/10/2024	PdxCre Tg/Azelia	Breeder
1L	F	10/28/2024	PdxCre +/tg	Breeder
1L1R	F	10/28/2024	PdxCre +/-	Breeder

Litters

#	DOB	Born	Survived	M/F	DoW
1	11/12/2024	3	3	2/1	01/15/2025

#	DOB	Born	Survived	M/F	DoW
2	01/08/2025	5	5	1/4	02/05/2025

Figure 2: Complete cage card with owner, protocol, animal table, litter table, and QR code

Phone QR workflow

1. Print the card at **100% scale**.
2. Use the **phone camera**, not a dedicated scanner.
3. Point at the **QR square**, not the 1D barcode.
4. Open the browser link.
5. Continue in either the workspace or the chat console.
6. Verify cage code and location before any write.

Common prompts used across all roles

The tutorial sections below show the exact seeded responses for these prompts:

- What needs attention today?
- Open cage F1-L01-C0006
- Show reports
- Print cage card for F1-L01-C0006
- Show project L01-PRJ-01
- What needs weaning this week?
- Show mortality follow-up
- Generate cage cards for Room A1
- Which labs are above expected load?
- What requests breached SLA?
- Show recent sample results
- Reserve 1 matching animal for project L01-PRJ-01




Figure 3: Phone scanning a printed QR cage card

- Generate a project closeout report

Role 1: Technician

What a technician cares about daily

A technician is trying to move quickly without losing biological or compliance accuracy.

Daily concerns usually are: - what cages need attention first - which cages are overdue for plug check or wean - whether counts and actual animals still match - whether a cage is blocked by protocol or welfare conditions - whether mortality, abnormal observations, or sample events were documented immediately

Typical technician data entries

Entry	Why it matters
male count / female count	reflects real room state and affects downstream planning
breeding status	changes expected tasks and litter timing
note text	captures abnormal findings at the moment they occur
task completion status	proves work was actually done
mortality record	supports welfare, veterinary, and census accuracy
litter / weaning values	affect survival metrics and future cage load
transfer destination	preserves chain-of-custody and location accuracy
sample collection state	affects genotyping and project readiness

Typical technician reports

- today's action list

- overdue task list
- cages with active alerts
- weaning list
- mortality follow-up list
- printable cage cards for the next room pass

Technician workflow in the workspace

1. Start from the dashboard or Cages

In the seeded dataset, the technician sees a room-facing operational surface with alerts, cage density, and due work.




Figure 4: Workspace view with cage alerts and room-facing operational context

Recommended first actions: 1. look for highlighted cages and alert badges 2. filter to the room you are entering 3. print or reprint cage cards if needed before rounds

2. Generate a cage card

Use the selected cage workflow to print F1-L01-C0006.

Expected printed card context: - Lab: Neurogenetics Lab - Location: Room A1 / Rack R1 - Protocol: IACUC-2026-0101 - Population: M3 / F3 / T6

3. Scan in the room and reopen the cage

After printing: 1. scan the QR 2. open the browser link 3. sign in if needed 4. continue from either mode




Figure 5: Workspace card-printing flow




Figure 6: Scan landing flow after QR open

4. Use the workspace for visual review

For F1-L01-C0006, the workspace helps the technician visually confirm: - strain and genotype summary - animal rows - litter timing - alert badges - whether a protocol or welfare issue is blocking edits

Technician workflow in chat

1. Start the shift with the morning brief

Prompt:

What needs attention today?

Observed response in the seeded training environment on **April 13, 2026**: - message: Here is the current operating brief. Start with the alerts and due tasks, then open a cage or report from here. Operational snapshot: - 6 active alerts - 6 high alerts - 6 tasks due now - 6 overdue tasks - 42 scoped cages - 84 scoped animals - top priority alerts included: - Protocol Expired on F1-L01-C0001 - Protocol Expired on F1-L01-C0002 - Protocol Expired on F1-L01-C0003 - top due tasks included: - Task #42 · wean on F1-L01-C0042 - Task #41 · plug_check on F1-L01-C0041 - Task #40 · wean on F1-L01-C0040

2. Open the same cage directly

Prompt:

Open cage F1-L01-C0006

Observed response: - message: Cage F1-L01-C0006 is open. You can inspect it, update counts, or add a note from here. Cage summary: - strain: Rosa26-LSL - genotype summary: fl/fl - population: M3 / F3 / T6 - location: Room A1 / Rack R1 - protocol: IACUC-2026-0101 - project: L01-PRJ-01 - DOB: 2026-03-06 - public scan page: /scan/tok_1_1_0006_4348

Tracked animals included: - F1-L01-C0006-SIRE - F1-L01-C0006-DAM - F1-L01-C0006-P01 - F1-L01-C0006-P02 - F1-L01-C0006-P03 - F1-L01-C0006-P04

3. Notice the abnormal condition before editing

This same cage had active alerts: - Protocol Expired - Open Vet Case - Protocol Deviation Open

That is an important technician lesson: **not every open cage is writable.**

4. Attempt the update and observe the hard-stop

Prompt:

Update cage F1-L01-C0006 males=3 females=4 status=Holding note=Tutorial walkthrough adjustment

Real result from the seeded dataset:

Protocol IACUC-2026-0101 expired on 2020-01-01

This is correct behavior. The app should stop room-side edits when protocol state blocks action.

5. Still do useful work: complete the due task

Prompt:

Complete task 6

Observed response: - message: Task #6 is now done. - updated task card: - task: #6 - type: wean - cage: F1-L01-C0006

6. Re-check the queue

Prompt:

Show overdue tasks

Observed response: - message: These are the current open tasks in your scope. Task pressure card: - 20 open tasks - 20 overdue - first rows in the open-task table included: - F1-L01-C0042 - F1-L01-C0041 - F1-L01-C0040 - F1-L01-C0039 - all four were still pending and assigned to Admin User

7. Ask for the weaning queue

Prompt:

What needs weaning this week?

Observed response in the refreshed seeded training environment on **April 24, 2026**:

- message: These litters are due or coming due for weaning in the next 7 days.

Weaning pressure:

- 11 litters due within 7 days
- 0 overdue

First rows in the weaning queue:

- F1-L01-C0015 · litter DOB 2026-04-03 · due to wean 2026-04-24 · 5 survived · Room A1 / Rack R1
- F1-L01-C0033 · litter DOB 2026-04-03 · due to wean 2026-04-24 · 5 survived · Room A1 / Rack R1
- F1-L01-C0012 · litter DOB 2026-04-04 · due to wean 2026-04-25 · 5 survived · Room A1 / Rack R1

8. Ask for mortality follow-up

Prompt:

Show mortality follow-up

Observed response:

- message: These mortality records still need necropsy or veterinary follow-up.

Mortality follow-up pressure:

- 20 open records
- 20 animals represented

First follow-up rows:

- #40 · cage F1-L02-C0037 · cause found dead · necropsy pending
- #39 · cage F1-L02-C0035 · cause found dead · necropsy pending
- #38 · cage F1-L02-C0033 · cause found dead · necropsy pending

The response also includes a direct Mortality CSV export link.

9. Batch-print the next room pass

Prompt:

Generate cage cards for Room A1

Observed response:

- message: Batch cage cards are ready for Room A1. The print view includes up to 40 cages so the URL remains browser-safe.

Batch print set:

- room: Room A1
- cards: 40

First included cards:

- F1-L01-C0001
- F1-L01-C0002
- F1-L01-C0003

The response includes a link like `/print/cards?ids=...` that opens the print view.

Technician overlap you should not skip

Even though these also matter to managers and researchers, a technician should practice them explicitly:

- printing the card - scanning with the phone camera - opening the cage from QR - checking alerts before editing - reading protocol state before changing counts - completing tasks in real time rather than later

Technician mini-mission

1. Print F1-L01-C0008.
2. Scan the printed card with a phone.
3. Open the cage in chat.
4. Ask `What needs weaning this week?`.
5. Ask `Show mortality follow-up`.
6. Ask `Show overdue tasks`.
7. Return to the workspace and visually confirm alerts for the same room.

Role 2: Facility Manager

What a facility manager cares about daily

A facility manager is less interested in one cage and more interested in whether the facility is controlled.

They care about: - room and rack utilization - which labs are above quota or trending up - protocol expirations and unresolved deviations - mortality and necropsy backlog - request SLAs and repeated breaches - which reports leadership, veterinarians, or IACUC will ask for next

Typical manager data entries

Entry	Why it matters
request status	keeps service flow and accountability visible
billing adjustments	affects chargeback accuracy

Entry	Why it matters
SLA settings	determines when delayed handoffs become actionable
protocol follow-up state	prevents non-compliant work from lingering
deviation / escalation state	supports corrective action and oversight
project priority / closeout state	helps resource allocation across labs
qualification or training updates	prevents blocked work from being invisible

Typical manager reports

- room utilization
- quota utilization by lab
- chargeback summary
- protocol expiration report
- mortality and necropsy report
- breeder productivity report
- survival report
- cohort handoff and closeout report

Facility manager workflow in the workspace

1. Start from analytics and compliance

The workspace is the strongest place to scan facility-wide signals quickly.




Figure 7: Compliance and operations visuals for facility-wide review

Use the workspace to review: - alert severity mix - protocol expiration watch list - deviation and compliance concentration - room and lab pressure at a glance

2. Move into analytics for capacity and flow

The seeded analytics workspace shows planning and throughput context that is hard to absorb in pure text.




Figure 8: Planner and analytics workspace used for manager review

Use this view to ask: - where is capacity tight - which cohort handoffs are stalling - which labs are breaching expectations repeatedly - what should be exported before the weekly meeting

3. Generate exports

From the workspace, a manager can generate: - closeout CSV/PDF - stalled handoff CSV/PDF - chargeback summaries - protocol usage and survival exports

Facility manager workflow in chat

1. Start with the morning brief

Prompt:

Give me the facility morning brief

Observed response in the seeded training environment on **April 13, 2026**: - message: Here is the facility morning brief. Alert load, due tasks, and scoped population are summarized first. Operational snapshot: - 6 active alerts - 6 high alerts - 6 tasks due now - 6 overdue tasks - 3,000 scoped cages - 372 scoped animals - priority alerts included multiple Protocol Expired cages in Neurogenetics Lab

2. Ask the next obvious question

Prompt:

Show protocol alerts

Observed response: - message: These protocols are expired or approaching expiration inside the alert window. - protocol alert table included: - IACUC-2026-014 · Synaptic Development Cohort · expires 2020-01-01 - IACUC-2026-0101 · Neurogenetics Lab Protocol 1 · expires 2020-01-01 - IACUC-2026-0201 · Synaptic Circuits Group Protocol 1 · expires 2020-01-01 - IACUC-2026-0202 · Synaptic Circuits Group Protocol 2 · expires 2020-01-01

3. Check room pressure

Prompt:

Show room utilization

Observed response: - message: Here is the current facility-level utilization, quota pressure, and chargeback snapshot. - room utilization table: - Room A1 occupied 3000 cages out of 240 capacity - utilization: 1250.0% - first rows in lab quota pressure: - Cancer Models Core at 89.6% - Aging Biology Lab at 119.87% - Renal Physiology Unit at 133.58% - Molecular Pathology Group at 133.09%

This is a good example of chat as a manager's quick exception console: concise numbers first, then the manager can switch to the workspace if they need a wider visual picture.

4. Pull chargeback or SLA context

Prompt:

Show chargeback summary

Observed response: - message: Here is the current facility-level utilization, quota pressure, and chargeback snapshot. - first rows in chargeback snapshot: - Cancer Models Core · 362 cages · 10,860 cage-days · estimated charge 9231.0 - Aging Biology Lab · 362 cages · 10,860 cage-days · estimated charge 9231.0 - Renal Physiology Unit · 362 cages · 10,860 cage-days · estimated charge 9231.0

Prompt:

Which labs are above expected load?

Observed response in the refreshed seeded training environment on **April 24, 2026**:

- message: These labs are above expected load or close enough to quota to deserve attention.

Load pressure:

- 5 labs above expected load
- 7 labs at or above 90%

First quota-watch rows:

- Renal Physiology Unit · 362 current cages · 271 expected · utilization 133.58%
- Molecular Pathology Group · 362 current cages · 272 expected · utilization 133.09%
- Pain Mechanisms Program · 152 current cages · 121 expected · utilization 125.62%

Prompt:

What requests breached SLA?

Observed response:

- message: Facility requests at or above 48 hours are treated as SLA breaches for this chat review.

Request SLA pressure in this refreshed seed:

- 0 open submitted/approved requests
- 0 breached SLA requests

Prompt:

Show reports

Observed response: - message: These are the reports and exports currently available from chat. - export links included: - Cages CSV - Cages Excel - Cages PDF - Breeder productivity CSV - Survival CSV - Protocol usage CSV - Mortality CSV - Cohort closeouts CSV/PDF - Stalled handoffs CSV/PDF - Billing statements CSV

Facility manager overlap you should not skip

Managers should still practice the same phone-based scan flow as technicians, because managers often verify conditions during rounds or incident follow-up.

Do not skip: - scanning a printed card - confirming a cage's protocol state on phone - using chat for a quick facility brief while away from a desk - returning to the workspace for a broader visual explanation

Facility manager mini-mission

1. In chat, ask Give me the facility morning brief.
2. Ask Show room utilization.
3. Ask Which labs are above expected load?.
4. Ask What requests breached SLA?.
5. Switch to the workspace and open analytics.
6. Export a cohort or handoff report.
7. Return to chat and ask Show protocol alerts.

Role 3: Researcher / PI

What a researcher cares about daily

A researcher or PI is usually asking: - do I have the right animals by genotype, sex, and timing - are cohorts being reserved and handed off correctly - are breeding lines supporting project demand - are sample and genotype results arriving fast enough - what story do the pedigree and closeout outcomes tell

Typical researcher data entries

Entry	Why it matters
project target counts	determines whether colony output meets study demand
genotype targeting rules	defines what animals are actually useful
cohort reservations / releases	links colony output to project planning
project priority or desired date	helps facility planning and breeder decisions
closeout notes and outcomes	captures what happened after animals were handed off
sample or genotyping state	ties line verification to experimental readiness

Typical researcher reports

- genotype-ready cohort list
- breeder productivity by line
- project cage list
- protocol usage report
- cohort closeout summary
- stalled handoff list
- sample and genotype result summary

Researcher workflow in the workspace

1. Use the workspace for samples, cohorts, and planning

The visual workspace is especially strong for researchers when they need to compare several things at once.




Figure 9: Samples and genotyping workspace with provider and result context

From the workspace, a researcher can: - review sample chain-of-custody - inspect provider workflow and results - review genotype-ready animals - apply or inspect project genotype target rules - reserve animals into a cohort - review handoff status and closeout history

2. Use pedigree when the biology matters

If a line behaves unexpectedly, open the pedigree explorer to see sire, dam, and pup relationships before making breeding or assignment decisions.

3. Use planner views when demand is changing

A researcher should understand not just today's available animals, but whether next month's demand is realistic.

Researcher workflow in chat

1. Open the main project

Prompt:

Show project L01-PRJ-01

Observed response: - message: Project L01-PRJ-01 is open. You can review handoffs, closeouts, and reporting from [here](#). Project card values: - title: Neurogenetics Lab Project 1 - lab: Neurogenetics Lab - target animals: 588 - assigned cages: 22 - active handoffs: 0 Handoff SLA card: - assigned max days: 2 - shipped max days: 5 - repeat breach threshold: 2 - source: default

2. Move from project to a specific cage

Prompt:

Open cage F1-L01-C0012

Observed response: - message: Cage F1-L01-C0012 is open. You can inspect it, update counts, or add a note from [here](#). Cage summary: - strain: C57BL/6J - genotype summary: Pending - population: M3 / F4 / T7 - location: Room A1 / Rack R1 - protocol: IACUC-2026-014 - project: L01-PRJ-02 - DOB: 2026-02-25 - tracked animals included F1-L01-C0012-P01 through P05 plus sire and dam - active cage alerts included Protocol Expired, Open Vet Case, Protocol Deviation Open, and Task Overdue

This seeded example is useful because it connects a project to a sample/result workflow.

3. Ask about reports or readiness

Prompt:

Show reports

Observed response: - message: These are the reports and exports currently available from chat. - report links included: - Cages CSV/Excel/PDF - Breeder productivity CSV - Survival CSV - Protocol usage CSV - Mortality CSV - Cohort closeouts CSV/PDF - Stalled handoffs CSV/PDF - Billing statements CSV

Prompt:

Show genotype-ready animals

Observed response: - message: These are the current genotype-ready animals and project matches in your scope. Genotype-ready snapshot: - 1 visible project - 20 ready animals in the returned table - 84 unassigned ready animals in scope Project readiness table: - L01-PRJ-01 in Neurogenetics Lab - 84 matching ready animals - 0 reserved - target 588 - recommended action: `assign_now` First ready animals returned: - F1-L01-C0042-P06 · f1/+ - F1-L01-C0042-P05 · WT/WT - F1-L01-C0042-P04 · tg/tg - each matched L01-PRJ-01 and was currently unassigned

Prompt:

Reserve 1 matching animal for project L01-PRJ-01

Observed response in the refreshed seeded training environment on **April 24, 2026**:

- message: Reserved 1 matching animal(s) for project L01-PRJ-01.

Reserved animal:

- F1-L01-C0042-P06
- sex F
- genotype fl/+
- cage F1-L01-C0042

This is a write action. It is intentionally explicit: the user must name the project and count.

Prompt:

Show stalled cohort handoffs

Observed response: - message: These are the current stalled cohort handoffs in your scope.
Stalled handoff snapshot: - 0 stalled assignments - 0 repeat-breach projects - 0 items in the 8d+ bucket -
the Stalled handoffs and Repeat breach watchlist tables were empty in this seeded run

4. Use the sample/result example

The tutorial-ready seed includes: - sample SMP-0004 - animal F1-L01-C0012-P01 - cage F1-L01-C0012

Use that chain to practice moving from colony context to sample context and back again.

Prompt:

Show recent sample results

Observed response:

- message: These are the most recent samples with received, resulted, or rejected states.

Sample result state:

- 5 recent review rows
- 2 resulted

First result rows:

- SMP-0011 · animal F1-L01-C0042-P01 · cage F1-L01-C0042 · status received · result/genotype WT/WT
- SMP-0008 · animal F1-L01-C0023-P01 · cage F1-L01-C0023 · status resulted · result/genotype fl/+
- SMP-0007 · animal F1-L01-C0017-P01 · cage F1-L01-C0017 · status received · result/genotype tg/tg

Prompt:

Generate a project closeout report

Observed response:

- message: Here is the closeout report shortcut. Use the links for exportable evidence, and the table for fast review.
- the refreshed seed had no closeout rows yet
- links included Cohort closeouts CSV and Cohort closeouts PDF

Researcher overlap you should not skip

Researchers should still be able to: - read a printed cage card - scan a QR from a phone during a lab conversation - switch from chat to the workspace when they need visual cohort or pedigree detail - understand why a protocol or welfare block changes the meaning of “available animals”

Researcher mini-mission

1. In chat, ask `Show project L01-PRJ-01`.
2. Open `F1-L01-C0012`.
3. In the workspace, inspect samples/genotyping.
4. Review `SMP-0004` and the linked animal context.
5. Return to chat and ask `Show recent sample results`.
6. Ask `Reserve 1 matching animal for project L01-PRJ-01`.
7. Ask `Generate a project closeout report`.

Shared Cross-Role Workflows

Workflow: Print a cage card

You can start this in either mode.

Workspace path

1. Open `Cages`.
2. Select a cage or filtered set of cages.
3. Choose `Generate + Print`.
4. Print at 100% scale.

Chat path

Prompt:

`Print cage card for F1-L01-C0006`

Observed response: - message: Cage card for F1-L01-C0006 is ready. Open the print view on a desktop or tablet, print at 100% scale, then scan the QR with a phone camera. Print

card returned: - badge: Breeding - lab: Neurogenetics Lab - location: Room A1 / Rack R1 - protocol: IACUC-2026-0101 - population: M3 / F3 Links: - print view: `/print/cards?ids=8` - scan page: `/scan/tok_1_1_0006_4348`

Workflow: Scan a printed card with a phone

1. Scan the QR square.
2. Open the browser link.
3. Verify cage code and location.
4. Continue in workspace or chat.
5. Save changes immediately.

Workflow: Move between modes without losing context

Recommended pattern: 1. use chat to open or update a cage quickly 2. jump into the workspace when you need lists, visuals, pedigree, or exports 3. come back to chat when you know the next direct action

This is not a compromise. It is the intended operating model.

Fun Applications To Learn Together

These scenarios make good team exercises because they combine biology, operations, and data quality.

1. Weekly colony huddle

Ask each learner to answer: - which cages are risky this week - which litters will drive next week's wean workload - which genotype results are still blocking the project

2. Facility operations check-in

Use the workspace to identify pressure visually, then ask chat for the concise version: - Give me the facility morning brief - Show stalled cohort handoffs - Show protocol alerts

3. Research planning conversation

Stand by a rack, scan a printed card, and ask: - is this cage actually useful for L01-PRJ-01 - what genotype do we still need - what breeder or sample action is next

4. Pedigree teaching moment

Open a breeding cage and compare: - what the card shows - what the pedigree view shows - what the genotype rule for the project requires

Troubleshooting

- Tutorial examples missing: you are likely using the scale-only seed instead of `seed_tutorial_demo.py`.
- Scan opens the wrong host: set `Scan Base URL` to a phone-reachable host or LAN IP, never `localhost` on printed cards.
- Phone camera cannot detect the code: make sure the printed card contains a QR square with enough contrast.
- Camera detects only a barcode app, not a browser link: scan the QR square, not the CODE128 bars.
- Report export fails: re-login and retry the endpoint.
- Card images missing: confirm the server is running and reload the print page.
- Edit blocked by protocol: the protocol is expired or otherwise not valid for direct cage edits.
- Same scan reopens repeatedly: refresh once after the latest build; scan tokens are now cleared from the URL after capture.

Final Checklist

A learner is ready to use Murisphere when they can: - read a cage card accurately - print and scan a cage card without help - explain when to use workspace versus chat - open a cage in both modes - recognize a protocol hard-stop - complete at least one role-specific write action - retrieve the report that matters for their role - explain why pedigree, genotype, litter timing, and protocol state all change operational decisions




Figure 10: Troubleshooting visual for QR scan and print quality